



Crystal Peeling & Cryo Modelling Excellence®

Crystal peeling is a comprehensive approach of skin rejuvenation that is non-surgical, safe and effective. Practises in Europe have been using dermabrasion for years. Now men and women are enjoying the unique results of this exciting procedure. With crystal peeling, you can see improvements in skin texture since the first session. It removes and peels off the dead cells of the epidermis as well as the pollution caused by daily exposure to dust. It corrects also all imperfections of the skin like stretch marks, acne scars, pigmentation and wrinkles. Oxynergy Paris device is one of the best device on the market to perform this treatment, medical CE certified and with 17 levels of abrasion.

- . A visible reduction of wrinkles and fine lines
- . A reduction of acne scars and stretch marks
- . A diminution of hyper-pigmented lesions (very efficient on hands)
- . An efficient treatment of photo ageing.
- . Results since the first session on the skin texture and softness
- . An improvement of collagen density and skin elasticity

Medical CE & Serial Number
2 years of warranty

© Oxynergy Paris 2008

WORKING PRINCIPLE:

This unique approach stimulates the production of fresh young skin cells, collagen and improves the Oxygen penetration into the skin. It is the ultimate advancement in non-surgical, non-invasive skin conditioning.

We perform a superficial skin polishing with the aid of a delicate hand piece, through which micro-crystals are projected. This technique delivers a gentle abrasion which efficiently produces satisfying results in as little as twenty minutes. Our device is also offering cryo therapy with an electronic cooling system to soothe the skin or administer skin care by temperature shock.

Features and advantages of EXCELLENCE CRYSTAL PEELING:

- Exclusive designed system
- Microprocessor controlled
- LCD negative blue screen
- Up to 5 m³ aspiration per hour
- Patented hand piece
- Disposable heads
- Low cost of consumibles
- Multi functional (Vacuum program and cryo-modelling)
- 17 levels of abrasion
- Auto cleaning system
- 2 years guarantee & CE medical
- Made in FRANCE under strict quality control
- Clinical study performed
- High quality abrasive powder (Al₂O₃)
- Peeling Hand piece
- Vacuum Massage Hand piece
- CRYO MODELLING Hand piece (-5°C to 15°C)

Important reminder on the efficiency of crystal peeling devices

- Highly purified quality of crystals particles
- Flexible power of the device
- Application system
- Potency of the air compressor