

Oxynergy Atomizer

The Oxynergy atomizer is an exclusive designed oxygen system concentrating all the advantages of the world newest technologies in Oxygen therapy. A unique and fashionable design, user-friendly operation. Oxygen therapy is providing sheer relaxation and impressive results on the quality of the skin. The Oxynergy Atomizer offers a wide range of applications to all skin types without any side effects. It's designed to be used with the full skin care range Oxynergy Paris.



- Oxygen **atomization**
- Oxygen **inhalation**
- **96% pure oxygen**
- Up to **6 liters per minute**
- 1 bar static pressure
- Electronic controlled system
- Integrated **oxygen quality level analyser**
- The **most silent oxygen generator** of the market < 30db
- **Made in germany** under strict quality control
- 2 years Warantee in factory
- Compact system (23 kg) easy to move (castor rolls)

CE - Serial Number - 2 years warantee

Working principle :

The atomizer draws in ambient air via a coarse dust filter and a suction filter, and passes the compressed air into containers with molecular sieves where the actual oxygen enrichment takes place. The oxygen contained in the air is then collected in the oxygen tank. and the outlet pressure regulated by a pressure reducer to 1 bar. Every time it is switched on, the appliance control system performs an internal self test of the entire system and the oxygen level. In operation the sequence control constantly monitors the electrical and pneumatic circuits of the system.

Components available:

- Flowmeter
- Inhalation masks & Oxygen humidifier
- Airbrush applicator
- Crystal sprayer
- Crystal needle

Indications:

- **Dermal application of skin care** with pure oxygen
- **Inhalation** of pure oxygen

Results:

- Improves skin microcirculation and oxygenation
- Reduces signs of skin ageing
- Improves freshness, tone and skin elasticity
- Reduces acne
- Skin soothing
- Relaxation & well being

