



# Oxygen Excellence®

OXYGEN EXCELLENCE by Oxynergy paris is an exclusive designed oxygen system concentrating all the advantages of the world newest technologies in Oxygen therapy. A unique and fashionable design, user-friendly operation.

OXYGEN EXCELLENCE offers a wide range of applications for all skin types without any side effects. With our complete skin care range you will achieve impressive results, provide wellness & sheer relaxation.

- . Skin and cellular oxygenation
- . Anti-ageing
- . Improvement of skin tone
- . Soothing effect on acne
- . A visible reduction of wrinkles and fine lines
- . Deep penetration of products
- . Long lasting results
- . High moisturizing treatment
- . Results since the first session on the skin texture and softness
- . An improvement of collagen density and skin elasticity
- . Improvement of skin microcirculation
- . Complementary to other medical treatments

CE & Serial Number  
2 years of warranty

© Oxynergy Paris 2008

## **WORKING PRINCIPLE:**

The device is producing pure oxygen drawing air through an exclusive system of 12 molecular sieves. Thanks to an electronic controlling system the oxygen is released at the desired continuous or pulsed flow. We are using hand made crystal handpieces to preserve the quality of the active ingredients and in order that they can be sterilized between each patient.

## **Features and advantages of OXYGEN EXCELLENCE :**

- Exclusive designed system
- Microprocessor controlled
- LCD negative blue screen
- Up to 15 litres pure oxygen per min.
- Crystal hand made sprayer and needle
- No consumibles
- Multi functional (4 Vacuum programs and 12 oxygen programs)
- 2 years guarantee & CE certified
- Made in FRANCE under strict quality control
- Clinical study performed
- Fast operation (no waiting time when changing program)

## **Treatment programs:**

### **ATOMIZATION**

Continuous flow of oxygen to spray serum over the skin, 4 adaptable levels 3, 6,10 and 15 litres.

### **MESO-OXYGEN**

Mesotherapy with oxygen to pulsed serums deep into the skin. The pulse frequency is 300 ms and 4 adaptable levels 3, 6,10 and 15 litres.

### **DIFFUSION**

Diffusion with oxygen to pulsed creams deep into the skin. The pulse frequency is 600 ms and 4 adaptable levels 3, 6,10 and 15 litres.

### **OXYENERGY MASSAGE**

Vacuum massage 4 programs of continuous and pulsed vacuum.

IMMEDIATE RESULTS  
WELLNESS & SHEER RELAXATION  
HIGH ADDED VALUED TREATMENTS  
COMPLETE PROFESSIONAL & SKIN CARE RANGE AVAILABLE  
FULL SUPPORT / MADE IN FRANCE / AFTER SALE SERVICE  
FULL SYNERGY WITH OTHER MEDICAL PROCEDURES